

Table S1. Primers for hybrid identification.

SSR	Forward primer (5'-3')	Reverse primer (5'-3')
<i>VVIP31</i>	TATCCAAGAGACAAATTCCCAC	TTCTCTTGTTTCCTGCAAATGG
<i>VMC8G6</i>	TCAGTAATCACGAGCTTCCCG	TGGAGTGGGGATATGGAAATG
<i>VRZAG83</i>	GGCGGAGGCGGTAGATGAGAGGGCG	ACGCAACGGCTAGTAAATACAACGG
<i>VVIN56</i>	GCCAAGTAGCCAAATTATAGAA	TTTATGCTCCGTGGTTTGAAAT
<i>VMC5H5</i>	TACTTGCCCAATGGGTAATGAC	ATGGCCTCCCTACAAAAGAAAC
<i>VVIF52</i>	AGGGAATTGAAGAGAAACTGTT	TCTGCCAAGCAAATGAAAGAAA
<i>VMC1G7</i>	AGCCATAAAGGCCTTAAAAAC	GGGTCCACATAGGTAGGAGATT
<i>VMC1G3.2</i>	ACTTAGCTTCAGAAGAAAATAGA	GATAGTTACCATACTTAGTCGGA
<i>VVIP31</i>	TATCCAAGAGACAAATTCCCAC	TTCTCTTGTTTCCTGCAAATGG
<i>VMC8G6</i>	TCAGTAATCACGAGCTTCCCG	TGGAGTGGGGATATGGAAATG