

**Supplementary TABLE 2. Adapters and primers used for MSAP analysis.**

<b>Adapters and primers</b>	<b><i>EcoR</i> I (E) (5'-3')</b>	<b><i>Hpa</i> II/<i>Msp</i> I (H/M) (5'-3')</b>
Adapter-1	CTCGTAGACTGCGTACC	GACGATGAGTCTCGAT
Adapter-2	AATTGGTACGCAGTCTAC	CGATCGAGACTCAT
Pre-amplification primers	GACTGCGTACCAATTCA	ATGAGTCTCGATCGGA
Selective amplification primers	GACTGCGTACCAATTCAAC(E32)	GATGAGTCTAGAACGGAGA(H39)
	GACTGCGTACCAATTCAAT(E34)	GATGAGTCTAGAACGGAGC(H40)
	GACTGCGTACCAATTCACA(E35)	GATGAGTCTAGAACGGAGG(H41)
	GACTGCGTACCAATTCAGT(E42)	GATGAGTCTAGAACGGATC(H44)
	GACTGCGTACCAATTCATG(E45)	GATGAGTCTAGAACGGCAT(H50)
	GACTGCGTACCAATTCATT(E46)	GATGAGTCTAGAACGGCTA(H59)
	GACTGCGTACCAATTCAG(E49)	GATGAGTCTAGAACGGCTC(H60)
	GACTGCGTACCAATTCGTC(E76)	GATGAGTCTAGAACGGCTG(H61)