

Supplementary Table S1 Primers used for qPCR assays

Gene	Forward primer (5' – 3')	Reverse primer (5' – 3')
<i>CchHSP10-01</i>	GAGACGGTGACTAAAGGTGGC	TCCAATCTTGACGCTTACAGG
<i>CchHSP20-02</i>	ACCGCCTTCATTTGCCGTTTG	TTCTCCCTCGCCTTCCCATCC
<i>CchHSP40-03</i>	AGGAAGGACTGAAAGGAGGAG	CTAAAGAATGAACCGAAGGGA
<i>CchHSP60-01</i>	GGTGCTATGAAGAAGATTGGT	AGGAGTTTATGAAGTACGGAGA
<i>CchHSP70-09</i>	AAGGAACGGGCAGATGCTAAA	TTGGTCCTCACCTCCTCATA
<i>CchHSP70-12</i>	TTCACCGATACTGAGCGACTC	ACGCAGGACATTCATACCAGA
<i>CchHSP90-02</i>	GTCTCCCGCATGAAGCCTAAC	AAACCACCACCTTCTCCACTT
<i>β-actin</i>	CTGGAAGGTGGACAGAGAGG	AAATCATCGCTCCACCAGAG