

Table S4 Primers used for transient dual-luciferase reporter

Gene	Forward Primer(5'-3')	Reverse Primer(5'-3')
<i>MiERF105</i>	CGGGATCCATGGCGTATCAAGA CGAAGCTTTTACC	GGGGTACCTCACATAACCGGGAGCTG AGG
<i>MiPAL2</i>	GGGGTACCGAATGAAATGAGC TTGACTAC	CGGGATCCTGCTTCACAATCAATCTCC CT
<i>MiANS3</i>	AACCTGGTTGGCTACTTTCTTCA GT	TGCCTCATTCAATAACAAATGAAGGGTAG G